GRAPHTEC

GS-CO2-UM-151

GS-CO2

For GL100

CO₂ Sensor

USER'S MANUAL

Thank you very much for buying this GRAPHTEC product.

This product is a measurement sensor (hereafter "module") that connects to the GL100-N/GL100-WL.

These directions describe preparations and cautions before measurement.

To ensure safety, please read the operation instructions, etc.

For details on the warnings and how to handle this module, please read the Quick Start Guide or USER'S MANUAL included on the CD-ROM (included in the GL100 packaging)

Confirmation of the exterior

After opening the package, please confirm that there are no problems (scratches and dirt) on the exterior before use.

Confirmation of the attached items.

User's manual (this book): 1

If by any chance faults are found, please contact the store where you bought the item.


* Please note that items mentioned in this book may change without prior notice.

604309071

MANUAL-CO

1 Part Names

This section describes the name and function of each part.



GL100 main module

	occurrence a main
2. Internal sensor	There is an internal sensor that detects the CO ₂ .
	The lamp for detecting flashes every two
	seconds.
3. Connector	Used to connect to the connector on the

GL100 module.

4. Cable packing . This packing is used when connecting the connector.

Used to mount to a wall

1 Hook portion

• The battery is not available for power supply. Use the USB cable to supply the power.

- Do not use the module close to your face. Your breath may
- affect measurements.

 Blocking the module's aperture may cause the airflow to

After connecting the GL100 to modules or sensors, please always check/set the

2 How To Measure

1. Power supply (Refer to Quick Start Guide or USER'S MANUAL.)

Supply power to this module via a USB cable connection.

WARNING This module cannot be powered with batteries.

Connecting the sensor will bring up the display "Sensor error!!" Remove the sensor, supply power to the module via USB cable connection and then operate it using the

2. Start-up and operation

(1) Screen display menu flow

After power-on, the GL100 is ready for operation by holding down [MENU] key. When the module is connected, "Module Type Recognition" screen is displayed. When the module is not connected, "Module Unconnected State" screen is

Operate in accordance with the displayed instructions.

FNSOR FRROR! "Please wait

Module unconnected state

<Operation> Connect the module

Recognition of module types



Press [ENTER] key

Module start-up

(2) Free-running screen



Hold down the [QUIT] key (approx. three seconds) to put the module into standby state. Press the [ENTER] key while in standby state to return to the free-running screen

3. Setting

(1) Setting screen operation

Item selecting screen

Press the [MENU] key on the free-running screen to go to the setting screen.

<How to set>

Select the item with the directional keys ($\triangle \nabla \triangleleft \triangleright$) and press the [ENTER] key.



If the submenu shows ↑ ↓ then there are selections in those directions

Numerical entry screen

<How to set>

Numbers can be inputted by increasing or decreasing the value with the \triangle and ∇ keys.



(2) DATA setting

Set the Sampling and Capture Mode those will be recorded to the data recording media.

The recorded data's size will be displayed in the information for the SD card being recorded to. Please take note of it.



DATA recording condition setting

Sampling	500 ms, 1, 2, 5, 10, 20, 30 s,	
	1, 2, 5, 10, 20, 30, 60 min	
Capture MODE	CONT, 1 Hour, 24 Hour	
Capture DIST	Memory, SD card	
	•	

(3) TRIGGER setting

Select the conditions for beginning data recording after measurement starts.

: Pressing the [START/STOP] key on this module will start/stop recording. : The recording will start with the trigger source conditions after

and will be stopped with the trigger source conditions.

pressing the [START/STOP] key. The recording will stop after pressing the [START/STOP] key. : The recording will start after pressing the [START/STOP] key



TRIGGER capture condition settings				
TRIG setting	Off, Start, Stop			
TRIG Source	Off			
	Alarm			
	Date	Date, Time		

(4) ALARM setting

The alarm information is required to set. Please set the number level when an alarm occurs.




ALARM settings				
Alarm	Off			
	Level	Off		
	/ Mode	Level		
		† H Value setting		
		↓ L * The level depends on the setting range.		

3 Recording


(1) Recording

Press the [START/STOP] key to start measuring with the set conditions.

After pressing [START] key, when the module is in awaiting recording start, "ARMED" is displayed, and then when recording is started, "REC" is displayed.



The module's status is shown with the lamp display.



CAUTION When accessing an SD card, do not remove the SD card. The data may not write properly or the SD card may be damag

(2) Recording completion

- Press the [START/STOP] key to stop measuring. • The screen display will change to the standby
- screen display.
- · Press [ENTER] key to change to the free-running screen display.



4 How To Confirm The Data

Check the recorded data with the application software included with this module using the method below (for details, refer to the USER'S

- (1) Connect the USB interface and check the online data
- (2) Insert the SD card into PC and check the data directly
- (3) Check the data directly from PC via wireless LAN

5 Specifications

Item	Contents		
Sensor method	NDIR system		
Measurement data	CO ₂ density		
Measurement channel	1 channel		
Measurement density range	0 to 9999 ppm		
Measurement density	±(5% of rdg + 30 ppm)		
accuracy	*Within the range of 0 to 5000 ppm		
Responsive 3 minutes or less			
Sampling interval	0.5, 1, 2, 5, 10, 20, 30 sec.		
	1, 2, 5, 10, 20, 30, 60 min.		
	* The processing of sensor is every 2 sec. Even if setting it to 0.5 or 1 sec, the data is processed at the same value (2 sec).		
Alarm	OFF / Level		
Cable length	approximate 20 cm		
Usage environment	0 to 50°C, 80% RH or less (non-condensing)		
External dimensions	66 x 100 x 24.6 mm (not including protruding parts)		
[W×D×H] (approximate)			
Weight (approximate)	99 g		


We recommend replacing the module periodically

6 Combined Measurement

1. Combined CO2 sensor and illumination / ultraviolet sensor measurement

Composite measurement can be done by using the branch adapter for GS (GS-DPA) and the illumination / ultraviolet sensor (GS-LXUV) (each sold

MARNING If the CO₂ sensor (GS-CO2) is included in the assembly, it cannot be powered with batteries.



(1) Screen display menu flow

After connecting the power supply, connect this module and operate it in accordance with the content displayed on the screen

- * Refer to "2 How To Measure" above and the USER'S MANUAL for the
- (2) Free-running screen




Hold down the [QUIT] key (approx. three seconds) to put the module into standby state. Press the [ENTER] key while in standby state to return to the free-running screen.

2. Combined CO2 sensor and temperature and humidity sensor measurement

Composite measurement can be done by using the branch adapter for GS (GS-DPA) and the temperature and humidity sensor (GS-TH) (each sold

MARNING If the CO₂ sensor (GS-CO2) is included in the assembly, it cannot be powered with batteries.




(1) Screen display menu flow

After connecting the power supply, connect this module and operate it in accordance with the content displayed on the screen.

* Refer to "2 How To Measure" above and the USER'S MANUAL for the temperature and humidity sensor.

(2) Free-running screen



Hold down the [QUIT] key (approx. three seconds) Press the [ENTER] key while in standby state to return to the free-running screen.

The module can be used approx. 1.5 m away from the GL100 by using an extension cable for GS (GS-EXC). However, you cannot connect and use multiple

It is not possible to be used by connecting two same sensors.

503-10 Shinano-cho, Totsuka-ku Yokohama 244-8503, Japan

GRAPHTEC August 1, 2014