



P11 - Rod



Full Specifications

PH Range:	0-14
Eo Point:	-25 to +25mV
Temp Range:	0-80°C
Response Time:	<30 secs, 95% FSD
Ref. Type:	AgCl, Liquid Filled, Single Junction
Junction Type:	Annular Ceramic
Shaft Length:	120mm \pm 2mm
Stem Diameter:	11.8 -12.2mm
Cable Length:	1 metre
Connector:	BNC (others available on request)
Combination:	Yes
Body Type:	Glass
Application:	Liquids, General Purpose with Toughened Membrane
Shelf Life:	12 Months

Accessories

Reference Refill Solution:	551-68 (150ml 3M KCl/AgCl)
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15