



### **Full Specifications**

PH Range:	0-14
Eo Point:	-25 to +25mV
Temp Range:	0-80°C
Response Time:	<30 secs, 95% FSD
Ref. Type:	AgCl, Double Junction, Liquid Filled
Junction Type:	Ceramic
Shaft Length:	120mm ±2mm
Stem Diameter:	11.8 -12.2mm
Cable Length:	1 metre
Connector:	BNC (others available on request)
Combination:	Yes
Body Type:	Glass
Application:	Liquids/Non Aqueous
Shelf Life:	12 Months

### **Accessories**

Reference Refill Solution:	551-55 (150ml LiCl)
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15