



Full Specifications

PH Range:	0-14
Eo Point:	0 to ± 25 mV
Temp Range:	0-100°C
Response Time:	<30 secs, 95% FSD
Ref. Type:	Silver/Silver Chloride
Junction Type:	Ceramic
Shaft Length:	120mm
Stem Diameter:	12mm
Cable Length:	1 metre
Connector:	BNC (others available on request)
Combination:	Yes
Body Type:	Glass
Application:	pH in photographic solutions and where sodium is present.
Shelf Life:	12 Months

Accessories

Storage Solution:	791-65
Recommended Cleaning Solution:	791-15