

## P11 Protected Bulb



### Full Specifications

PH Range:	0-14
Eo Point:	0 to $\pm 25$ mV
Temp Range:	0-100°C
Response Time:	<30 secs, 95% FSD
Ref. Type:	Ag/AgCl, Single Junction
Junction Type:	Porous Teflon
Shaft Length:	120mm
Stem Diameter:	12mm
Cable Length:	1 metre
Connector:	BNC (others available on request)
Combination:	Yes
Body Type:	Glass
Application:	pH measurements where a very robust device is required
Shelf Life:	12 Months

### Accessories

Reference Refill Solution:	551-63 (3.5M KCl/AgCl Partial Gel)
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15